

Product Overview

The **DeNovix DP-02 Pipette** features a comfortable, ergonomic design and delivers 0.1 – 2.0 μL volume aliquots. The pipette is a perfect companion for the DS-11 Spectrophotometer.


Each pipette may be fully disassembled and can be autoclaved. The DP-02's adaptable design allows it to accept microvolume tips from all major suppliers.



Quality Assurance

Each pipette is assembled, calibrated and certified according to EN ISO 8655. Each device is supplied with a quality control certificate from our ISO 9001 registered manufacturing facility.

Features



Specifications

Volume (μL)	Accuracy (%)	Precision (%)
0.2	± 12.0	± 12.0
1.0	± 2.7	± 1.3
2.0	± 1.5	± 0.7