

Squid[™]

Full Range Pipette (1 - 1000 μ L)



DeNovix[®]

A Full Set of Pipettes in One Device: 1 - 1000 μL

Squid™ Pipette: Simplify Your Bench

- Replaces a Full Set of Pipettes
- 1 - 1000 μL in a Single Device
- Universal Tip Compatibility
- One-Click Switching Between Tip Sizes

How Does It Work?

The unrivaled performance of Squid is achieved through **Dynamic Volume Control™**, which combines a patented tip selector mechanism and DeNovix expertise in microvolume motor control.

Squid employs two tip sizes, spanning the full 1 - 1000 μL range to eliminate the need to switch devices mid-workflow. Instantly activate each mode with a single press of the trigger button.

Switch with the Press of a Button



Large Tip Mode
100 - 1000 μL



Small Tip Mode
1 - 200 μL



Designed by Life Scientists

- **Ergonomic Design:** A comfortable, balanced grip for prolonged use
- **Custom Presets:** Instantly switch between required volumes
- **Lab-Ready:** Pre-configured modes and custom options available
- **All-Day Battery Life:** Equipped to support your workflow

Exceeds ISO Calibration Specifications Across the Full Range

The Squid Pipette outperforms industry standards across its entire 1 - 1000 μL range, doing the work of four traditional pipettes without compromising on data quality.

Squid Performance Against ISO 8655		
ISO Range	Systematic Error (Accuracy)	Random Error (Precision)
1000 - 100 μL	✓	✓
200 - 20 μL	✓	✓
20 - 2 μL	✓	✓
10 - 1 μL	✓	✓



View Performance Data

Tested and certified to outperform ISO 8655 requirements from 1 μL to 1000 μL .

A Full Set in One Pipette

The **DeNovix Squid™** is a revolutionary electronic pipette that replaces an entire traditional set—from **P10** to **P1000**—with a single, high-performance device.



Personalize Your Lab with 5 Color Choices



Accessories

DeNovix Universal Pipette Tips

20 μ L Pipette Tip	Neutral, Sterile, Filtered, Racked
200 μ L Pipette Tip	Neutral, Sterile, Filtered, Racked
1000 μ L Pipette Tip	Neutral, Sterile, Filtered, Racked
200 μ L Pipette Tip	Neutral, Sterile, Racked
1000 μ L Pipette Tip	Neutral, Sterile, Racked
200 μ L Pipette Tip	Neutral, Racked
1000 μ L Pipette Tip	Neutral, Racked

Squid Charging Stand

Charging Stand SKU: A-70300

Calibration and Maintenance Services

The Squid Pipette is fully compatible with standard gravimetric calibration procedures performed by any accredited pipette calibration provider, integrating seamlessly into existing laboratory calibration schedules. For customers in the USA, DeNovix also offers a calibration service to provide a convenient option backed by the expertise of the instrument's manufacturer.

Specifications

Pipetting Type	Air Displacement Pipette
Volume Range	1 - 1000 μ L
Accuracy / Precision	Exceeds ISO 8655-2_2022 at All Volumes
Warranty	2 Years
Weight	230 grams
Autoclavable	Yes - Tip Holders and Ejectors
Dimensions (L x W)	30.4 x 4.4 cm
Display	1.5 Inch Full-Color HD Screen
Connectivity	Bluetooth 5.0, Wi-Fi 2.4 GHz
Power Supply Input	100 - 240 V ~ 50/60 Hz
Power Supply Output	5 V / 3.0 A / 15.0 W
Battery Type	1200 mAh Li Ion UL-1642 Compliant
Charging	USB-C Cable (Included) Squid Charging Stand (Optional)
Approvals	CE, UL/CSA, FCC and Japan CAB/VCCI
Patents	denovix.com/patents



+1.302.442.6911 | info@denovix.com
Copyright ©2026 DeNovix Inc. | 2-APR-2026

Order Now